



BAG FILTERS

Item No.	Description
GA-600	Small Fountain Solution Recirculation System Filters
GA-602	Large Open Mouth Fountain Solution Recirculation System Filters

Packaging: Per Piece

Application: These bag filters are a replacement for more cumbersome canister filters. Some of the canister types contain a short-lived foam filter. The bag is double-lined and contains an extra filter media inside.

Advantages:
Effectively Filters Emulsified Inks
Keeps Contaminants From Circulating
Very Easy to Install

Usage & Recommendations:

How to use: Installation is quick and easy. Disconnect your canister or other filter from the tank's incoming water line. Slide the incoming water line over the ½ inch barb on the end of the Gans filter. (See white barb in photo.)

When to replace: When your fountain solution water begins to look dirty, just pull the used bag filter off and slide a new one on.

To Optimize Performance, Use These Related Products:

GA-503, Gans System Flush

